



*This simple to use kit delivers a professional solution featuring a convenient, pre-portioned disinfecting cleaner on the EPA List "N".*

**A Professional Solution:**

- Tested and confirmed 30-day efficacy once mixed and added to the wipes
- Features an EPA-approved, "N" listed disinfectant cleaner
- 5 minute kill claim for pathogens (emerging or otherwise)

**Features:**

- Can be used to clean and disinfect all hard surfaces, equipment, tools and more.
- Wipers easily absorb the cleaning solution while remaining durable to handle your daily cleaning needs.
- Reusable bucket has a *custom designed, airtight, flip top lid* to keep your wipes saturated, avoid contamination and reduce staff exposure to the chemicals.
- Each coreless roll pulls from the center – perfect for feeding wipers through the kit's bucket lid.

**(800) 615-8699**

**www.liquidsafety.com**



Item#	Description	Count
-------	-------------	-------

94400EF	Disinfecting Wipe Kit Pro: Wiper Roll (300 Wipes) & Disinfecting Solution	Each
74491EF	Disinfecting Wipe Kit Pro Refill: Wiper Roll & Disinfecting Solution	6 Refills/Case




**Industrial**



**Retail**



**Medical**



**Restrooms**

**Included with each kit and refill:**

***128 E-FECTicide Disinfecting Cleaner***

- One-step cleaning and disinfecting solution
- EPA-approved for Emerging Pathogens (List "N" EPA Reg No. 6836-365-5449)
- **Convenient, pre-portioned, easy-to-mix concentrate**, each 2 oz bottle makes enough for one kit (64 OZ)
- **Professional** quat-based solution tested in the wipe kit for 30 days

Our Chemical Partner:

**Multi-Clean®**

*Now with a custom designed, airtight, flip top lid!*



## Technical Specs:

### Disinfecting Wipe Kit Pro

- FDA (HBD) Approved White Bucket w/Lid (7.375" High x 7.8125" in Diameter)
- Dry Wiper Roll: 300 Airlaid Cellulose-Based Wipers (6.8" x 6")
- Disinfecting Concentrate
- Instruction/OSHA Sheet & Bucket Label
- 150 Kits Per Pallet



### Disinfecting Wipe Kit Pro Refill

- 6 Dry Wiper Rolls: 300 Airlaid Cellulose-Based Wipers (6.8" x 6") Per Roll
- 6 Bottles of Disinfecting Concentrate
- Instruction/OSHA Sheet
- Shipping Box: 13.3" x 20" x 7.6"
- 84 Cases (of 6 rolls) Per Pallet



### Emerging Pathogens Claim:

*Multi-Clean 128 E-FECTicide Disinfectant Cleaner (EPA registration No. 6836-365-5449) has demonstrated effectiveness against viruses similar to the 2019 Novel Coronavirus (SARS-CoV-2, the cause of COVID-19) on hard, non-porous surfaces. Therefore, Multi-Clean 128 E-FECTicide Disinfectant Cleaner can be used against the 2019 Novel Coronavirus (SARS-CoV-2, the cause of COVID-19) when used in accordance with directions for use against Norovirus on hard, non-porous surfaces.*

*This claim is based on compliance with product instructions and a 5-minute wet time.*

Our Chemical Partner:

**Multi-Clean®**

### Ideal For:

Manufacturing Environments	Educational Facilities & Campuses
Commercial Buildings	Retail Stores
Health Care Settings	Day Cares
Hospitality	Contract Services
Property Management	Food Service Operations

## Quick Guide

**Dilution** 1:64 (1370 ppm quat minimum)

One 2 oz bottle makes enough for one Disinfecting Wipe Kit Pro (64 oz). Pour the entire mixture over the wipers (in the bucket) in a circular motion from the outside in to the center.

### Uses:

For cleaning and disinfecting non-porous high touch surfaces. For heavily soiled areas, a pre-cleaning step is required. Professional grade disinfectant approved for use to disinfect against SAR-CoV-2, the cause of COVID-19. EPA Reg No. 6836-365-5449

### Surfaces:

Use on any non-porous surface not harmed by water. Countertops, ceramic, painted metal, stainless steel, plastic, naugahyde, sealed wood.

### Restrictions:

Use good safety practices when handling concentrated disinfecting solution. *See SDS for further information.*

 **Spilfyter® Disinfecting Wipe Kit Pro**

**LIQUID SAFETY SOLUTIONS**  
BY FYTERTECH NONWOVENS